




BMC iVolve Full Face Mask Sizing Guide

Choosing the correct cushion size

- There are three different mask sizes available. Small (S), Medium (M) and Large (L)
- The length (L1) is measured as shown below.
- Choose the smallest size mask cushion that is long enough.



BMC	
Full Face Mask Measurements	
(S) SMALL	= 8 ~ 9cm
(M) MEDIUM	= 9 ~ 10cm
(L) LARGE	= 10 ~ 11.5cm



IMPORTANT NOTES

- For the best results, we highly recommend a personal fitting.
- We understand that this is not always possible, so we have provided this information to assist.
- For online purchases and phone orders, we do not accept responsibility if the mask is uncomfortable or does not fit.
- For hygiene reasons, once the sealed mask packet has been opened we cannot exchange it.
- For online purchases and for phone orders, you agree to accept the limitations described above.
- If you do not agree to these limitations, please do not purchase.