




BMC Micro Nasal Mask Sizing Guide

Choosing the correct cushion size

- There are three different mask sizes available. Small (S), Medium (M) and Large (L)
- The length of the nose is “N1”, and the width of the nose is “N2”.
- Choose the smallest size mask cushion that is wide enough to extend beyond the nostrils but does not obstruct breathing.



For the different sizes, please refer to table below.

Dimensions of Face and the Mask

		Small = S (mm)	Medium = M (mm)	Large = L (mm)
BMC Micro Nasal Mask	Length (N1)	52	56	61
	Width (N2)	48	51	53

IMPORTANT NOTES

- For the best results, we highly recommend a personal fitting.
- We understand that this is not always possible, so we have provided this information to assist.
- For online purchases and phone orders, we do not accept responsibility if the mask is uncomfortable or does not fit.
- For hygiene reasons, once the sealed mask packet has been opened we cannot exchange it.
- For online purchases and for phone orders, you agree to accept the limitations described above.
- If you do not agree to these limitations, please do not purchase.